

bestDUAL®

LESS INVASIVE SPECIMEN ACQUIREMENT SYSTEM

BONE MARROW TREPINE BIOPSY NEEDLE WITH **RETENTION TRAP** AND **HANDLE GRIP ADAPTING TO HAND SIZE**



HANDLE GRIP ADAPTING TO OPERATOR'S HAND SIZE

NOLOSE® SYSTEM

The inner sample retaining system that always guarantees success in specimen acquirement, free of damage.

NOLOSE® SYSTEM

SAMPLE RETENTION TRAP



- Cylindrical expansion chamber
- Cylindrical entrance with TREPAN® sharpening
- High penetration QUADRAN® bevel

COMPARING SPECIMEN'S DIAMETER

8G Jamshidi = Ø 2,4 mm | 9G BEST LISAS DUAL = Ø 2,5 mm
11G Jamshidi = Ø 1,7 mm | 11G BEST LISAS DUAL = Ø 2,1 mm

BONE MARROW TREPINE BIOPSY NEEDLE WITH RETENTION TRAP AND HANDLE GRIP ADAPTING TO HAND SIZE

Order Codes:

9G x 10 - 15cm = BSL0910D

11G x 10 - 15cm = BSL1110D

13G x 5 - 9cm = BSL1309D

The advanced design of BEST-LISAS DUAL® will guarantee the following advantages:

- ▶ **Variable ergonomic grip handle**
BEST-LISAS DUAL® has a variable grip handle, suited to a **large (man's) hand** or a **small (woman's) hand**. The grip handle is also fitted with Luer Lock.
- ▶ **Ease of use and guaranteed specimen collection**
The unique NOLOSE® core retention system assures success for a majority of cases.
For soft bones or critical pathologies the use of BEST-LISAS DUAL® needle combined with the accessory SAFE-LOCK® will always guarantee success in the specimen acquisition.
NO MORE LOST CORES!
- ▶ **A considerable reduction in pain**
Without any Jamshidi-type point (conical tapered), BEST.LISAS DUAL is considerably less traumatic.
- ▶ **A more compact and abundant specimen**
The TREPAN® crown edge with its advanced geometry will guarantee an intact and larger specimen core even in the softest bone with less pain to the patient and ease of use for the physician.
- ▶ Colour coded for easier needle size (gauge) determination.
- ▶ Sterile, single use: reduction of infection hazard.

bestDUAL®

AVAILABLE SIZES

SIZE G (mm)	LENGTH cm	COLOUR CODE	PCKG pz-units
9 (3,50)	10 - 15	pink	20
11 (3,00)	10 - 15	white	20
13 (2,40)	5 - 9	light blue	20



Address: 31 Airlie Road, Pullenvale, Queensland 4069 Australia
Phone: 07 31669769 Fax: 07 30111029 Mobile: 0410 537 746
Email: sales@newmedical.com.au Web: www.newmedical.com.au
ABN: 63 124 888 348

